



North America


Infinity-Lock™ Button System

For fixation during the healing process of acromioclavicular separations due to coracoclavicular ligament disruption.


Safe
Strong
Simple

Why Infinity-Lock™ Button System?

- Use the Infinity-Lock Button System for AC joint separation requiring a simple and replicable operative solution.
- Woven polyester implant utilizes medical manufacturing techniques supported by over 30 years of data.
- Device and surgical technique designed by leaders in upper limb surgery.
- One size fits all, eliminating potential for mis-sizing and stocking issues.
- Titanium alloy Button provides high strength, low profile fixation on the clavicle.
- Infinity-Lock Button System is a strong and reliable alternative ACJ repair method.



How does it work?



The information above MUST NOT be used to perform the surgical procedure.

Always refer to the Instructions for Use LAB 272 which contains detailed information and essential warnings and precautions.

Order codes

102-1089 Infinity-Lock Button System (supplied sterile) includes:

- Infinity-Lock Tube-Tape, 7 mm x 240 mm
- Infinity-Lock Button, 4 mm x 12 mm
- Cannulated drill bit, plain shank to fit Jacobs chuck, 4.0 mm diameter x 120 mm
- Guidewire, diameter 2.0 mm x 150 mm

Optional Disposable Instrument (supplied sterile):

202-1411 CC-Hook, with a curved end, Left

202-1413 CC-Hook, with a curved end, Right

Contact **Xiros™** Inc. on enquiries@xirosna.com for further details or visit our website at www.xirosna.com

We are registered with the FDA in the USA, Registration number 3015525389, FEI Number 3015525389.

Developed and manufactured by **Xiros™** Ltd.

Springfield House, Whitehouse Lane, Leeds LS19 7UE

Xiros, and Infinity-Lock are trademarks of Xiros. All rights reserved © Xiros™ 2024. Worldwide patents and patents pending.

Get involved

At Xiros we are constantly striving to improve outcomes for patients, and we want to hear from you.

Help us create the future of sports medicine and share your opinions and expertise with us at:

www.xirosna.com